

2x AbTaq Master Mix

Catalogue No.	Description	Pack Size
OP-MM-001-100	100 reactions of 20 µL	1 x 1 mL
OP-MM-001-500	500 reactions of 20 µL	5 x 1 mL

Product Description

2x AbTaq Master Mix is an **efficient, fast and convenient** way to perform routine PCR. Our engineered AbTaq polymerase supports robust PCR amplification over wide range of samples and even in the presence of commonly known PCR inhibitors. The ready-to-use 2x Master Mix is designed to minimize handling, yet provide optimal AbTaq polymerase activity and PCR performance. The generated PCR product can be directly loaded into the gel without addition of DNA loading buffer.

Product Application

The 2x AbTaq Master Mix is suitable for PCR amplification from the following templates:

- cDNA
- purified plasmid, genomic DNA
- crude/impure DNA sample preparations

and is ideal for

- PCR-based fingerprinting including but not limited to genotyping, colony PCR etc.

Features	Advantages
AbTaq 2.0 Polymerase	<ul style="list-style-type: none"> • Amplification of up to 4 kbp long fragment from a plasmid template • Fast and reliable extension speed of 15 sec / kbp Resistant to high concentrations of common PCR inhibitors e.g.: <ul style="list-style-type: none"> - up to 1 mM EDTA - 200 mM urea - >10% fetal bovine serum - up to 5% EtOH thus reducing the need to obtain cleaner DNA samples
2x Master Mix	The 2x concentrate contains all the components required for a PCR reaction, except water, DNA template and primers. It is formulated to enable (i) optimal amplification with our AbTaq polymerase and (ii) direct loading on the gel. With less pipetting steps required for PCR reaction setup and analysis, the 2x

	AbTaq Master Mix saves time, materials and increases assay reproducibility.
--	---

Product Protocol

PCR reaction

The following PCR reaction conditions are suggested as a starting point for further optimization.

Component	Final concentration	PCR reaction-volume	
		20 µL	50 µL
2x Master Mix	1x	10 µL	25 µL
Forward primer (10 µM)	400 nM	0.8 µL	2 µL
Reverse primer (10 µM)	400 nM	0.8 µL	2 µL
Template	variable	variable	variable
ddH ₂ O		to 20 µL	to 50 µL

PCR amplification conditions

Step in PCR cycle	Temp. °C	Duration	Cycle(s)
Step 1 : Initial denaturation	95	2 min	1
Step 2: Denaturation	95	15 sec	} 33
Step 3: Annealing	50-65	15-30 sec	
Step 4: Extension	72	15 sec / kbp	
Step 5: Final extension	72	2 min	1
Step 6: Hold	10	∞	1

Important considerations

PCR Reaction

Our suggested reaction **volume** of 20 µL is sufficient for the vast majority of assays but can be scaled up or down if necessary. For multiplex PCR, we do not recommend using < 20 µL reaction volume.

Primers should be designed according to standard PCR guidelines and are typically used in the concentration range of 100-400 nM.

For the amount of **template** to use, a wide range, from 0.1 pg to 1 µg can be used successfully for a PCR amplification. The optimal amount will vary depending on the quality and complexity of the DNA sample preparation. In a 50 µL PCR reaction, our recommended amount is (i) 10 pg - 50 ng for plasmids, and viral genomes, (ii) 10 ng - 200 ng for complex prokaryotic and eukaryotic genomes, and (iii) 2 µL from a standard 20 µL cDNA reverse transcription reaction converting up to 2 µg of RNA.

For problematic **templates with high GC** content, up to 10% of DMSO is recommend to be added to a PCR reaction to relax the DNA secondary structures.

PCR Amplification Conditions

Initial denaturation of the DNA strands is performed at 95°C for 1 min for simple DNA templates (cDNA, plasmid) and can be increased up to 3 min for complex eukaryotic genomes or templates rich in GC. For colony PCR assays, 2 - 4 min is suggested to ensure complete bacteria lysis. For cycling, we recommend that DNA **denaturation** is set at 95°C for 15 s.

The **annealing temperature** will depend on the user's primer design and is usually 3 - 5 °C below the lowest primer melting temperature (T_m). Ideally, the optimal annealing temperature should be empirically determined by performing a temperature gradient PCR in the range +/- 5 °C of the calculated T_m of the primers. If DMSO is added to the PCR reaction, the annealing temperature should be reduced by 0.75°C for every 1% of DMSO added (by volume).

The **extension** step is performed at 72°C and the duration will be dependent on the length of the PCR product and the template complexity. For low complexity templates (e.g. cDNA, plasmid), the recommended extension time is 1 kbp / 15s. The duration should at least be doubled for genomic DNA.

Typically, a PCR reaction will consist of 25-35 **cycles**. Up to 40 cycles may be necessary to enable sufficient amplification of a low-abundance template.

The **final extension** step is set to fill any incomplete ends of the newly synthesized strands and is performed at 72°C. Duration of this step will dependent on the length and composition of the PCR fragment but 2 - 10 min is usually sufficient.

Storage Condition and Stability

Stored at -20°C. Avoid repeated freezing / thawing cycles. For handling smaller sample sizes, we recommend aliquoting the 2x Master Mix before storage.

Technical support

If you are experiencing difficulties with setting up your PCR reaction and/or PCR amplification conditions, please contact our team with any relevant information at infoab@abasiabiolabs.com