

# AbScript 3.0 Reverse Transcriptase

Catalogue No.	Pack Size	Concentration
SP-RT-001-2000U	2000U	200U/ul
SP-RT-001-10K	10000U	200U/ul
SP-RT-001-40K	4 x 10000U	200U/ul

## Product Description

AbScript™ 3.0 Reverse Transcriptase is an enhanced M-MLV RNA-dependent DNA polymerase capable of high-quality long cDNA synthesis (> 9kDa), increased thermostability and reduced RNase activity. Synthesis of cDNA is possible at temperatures up to 55°C.

## Product Contents

Catalog No.	Component Name	Size
SP-RT-001-2000U	AbScript 3.0 Reverse Transcriptase	2000U
	5x AbScript Buffer	1 x 1ml
	0.1M DTT	1 x 0.5ml
SP-RT-001-10K	AbScript 3.0 Reverse Transcriptase	10000U
	5x AbScript Buffer	1 x 1ml
	0.1M DTT	1 x 0.5ml
SP-RT-001-40K	AbScript 3.0 Reverse Transcriptase	4 X 10000U
	5x AbScript Buffer	4 x 1ml
	0.1M DTT	4 x 0.5ml

## Storage Condition and Stability

The enzyme and undiluted buffer solutions are stable when stored at -20°C for 1 year.

## Unit Definition

One unit is defined as the amount of enzyme that will incorporate 1 nmol of dTTP in 10 minutes at 37°C, using poly(A) and oligo(dT) as the template and primer, respectively.

## Product Application

RT-PCR, cDNA synthesis, qRT-PCR

Features	Specification
Concentration	200U/ul
Recombinant Enzyme	Yes, an RNA dependent DNA polymerase
Enzyme Reaction	Reverse Transcription; first strand synthesis
Sample/Target	RNA

Included dNTPs	No
Reaction Condition	At 1x concentration the reaction buffer assures optimal activity for AbScript 3.0 Reverse transcriptase

## Product Suggested Procedure

1. Thaw the reagents and store on ice. Briefly vortex and centrifuge all reagents before setting up the reactions.
2. In a sterile reaction tube on ice, add the following components listed below:

Components	Volume	Concentration
<b>Primer</b>	1uL	50-250ng of random primers; or 2 pmol of gene-specific primers
<b>RNA</b>	X uL	10pg- 5ug of total RNA or 10pg – 500ng of mRNA
<b>10mM PCR Grade Nucleotide Mix</b>	1uL	-
<b>Water, PCR Grade</b>	To make a final vol of 13ul	-

3. Mix and centrifuge briefly. Heat 65°C for 5 minutes. *Optional: If using random primers, incubate tube at 25°C for 5 minutes before proceeding with incubation at 50°C for 30-60 minutes. Additionally, increase the reaction temperature to 55°C for gene specific primers.*
4. Incubate on ice for at least 1 minute.
5. After incubation, add the following components listed below to a **final volume of 20ul**:

Components	Volume	Working concentration
<b>5x AbScript Buffer</b>	4uL	1x
<b>0.1M DTT</b>	1 uL	500 uM
<b>RNase Inhibitor</b>	1 uL	40U
<b>AbScript 3.0 Reverse Transcriptase</b>	1 ul	200U. Adjust up to 2uL (400U) if increased yield is needed.

6. Incubate at 50°C for 30-60 minutes.
7. Denature AbScript 2.0 Reverse Transcriptase by heating at 70°C for 15 minutes.