

AbTaq™ 2.0 Polymerase

Catalogue No.	Pack Size	Concentration
OP-TP-001-200U	200U	5U/ul
OP-TP-001-1000U	1000U	5U/ul
OP-TP-001-2500U	2500U	5U/ul

Product Description

AbTaq™ 2.0 Polymerase is a Taq DNA polymerase enzyme engineered for fast extension speeds, thermal stability and a hot start function. AbTaq™ 2.0 Polymerase is tolerant of common PCR inhibitors such as whole blood.

Product Contents

Catalog No.	Component Name	Size
OP-TP-001-200U	AbTaq 2.0 Polymerase	200U
	10x AbTaq Buffer	1 x 1ml
OP-TP-001-1000U	AbTaq 2.0 Polymerase	1000U
	10x AbTaq Buffer	1 x 5ml
OP-TP-001-2500U	AbTaq 2.0 Polymerase	2500
	10x AbTaq Buffer	2 x 7ml

Storage Condition and Stability

The enzyme and undiluted buffer solutions are stable when stored at -20°C for 1 year.

Unit Definition

One Unit will incorporate 10 nmole of deoxyribonucleotide into DNA in 30 minutes at 75°C.

Product Application

RT-PCR, genotyping, PCR, qRT-PCR

Specification	Features
Concentration	5U/ul
Recombinant Enzyme	Yes
Enzyme Activity	5' →3' exonuclease activity
Sample/Target	Genomic DNA, plasmid DNA, cDNA etc.
With Hotstart	Yes
dNTPs included	No
Reaction Condition	At 1x concentration the reaction buffer assures optimal activity for AbTaq™ 2.0 Polymerase

Product Suggested Protocol

1. Thaw the buffer and store on ice. Briefly vortex and centrifuge all reagents before setting up the reactions.

For research use only (RUO).
Not for use in diagnostic procedures

2. In a sterile PCR reaction tube placed on ice, add the following components listed below:

Components	Volume	Final Concentration / Notes
Primer Fwd	Variable	200 nM
Primer Rev	Variable	200 nM
10x AbTaq reaction buffer	2.5ul	1x
10mM dNTP mix	0.5ul	200 µM
DNA template	Variable	0.1 – 100ng
2x AbTaq Accelerator Cocktail (optional)	12.5ul	Sold separately. May be replaced with water if not required
AbTaq 2.0 Polymerase at 5U/ul	0.1ul	0.5 – 2 U, adjusted based on user optimisation
Water, PCR Grade	Variable	To a final vol of 25ul

3. Mix and centrifuge briefly.
4. Suggested cycling conditions for short PCR products (<600bp):

Duration	Temp	No. of cycles	Optimisation Notes
1 min	95°C	1 cycle	Time may be extended to 8 min if crude samples containing blood/plasma/serum are used
25 sec	95°C	25 – 40 cycles	40 – 60 sec
45 sec	50 – 68°C		40 – 60 sec, temperature to depend on primer Tm
24 sec	68°C		To be adjusted based on an extension rate of 2kb/min

5. Run and visualise PCR products on a suitably stained agarose gel.